

# Safe Biopsy Container

Art.Code: 1748, 1749 Document type: Sales specification

CE IVD

ISSUED BY: Raphaël Roudet

DATE: 2023-02-25

DOC ID

REVISED BY: -

DATE: N/A

202302251520

APPROVED BY: Jonny Jansson

STATUS: APPROVED

## THE ROAD TO SAFE BIOPSY HANDLING.

**Eliminate your staff exposure to formalin with Biopsafe®, an innovative solution for safe transport of biopsy samples between department and their storage.**

The smart system consists of a practical small and easy to use container with formalin encapsulated in the lid. The release of the formalin is done hermetically with a press of the thumb and the safety of medical staff is guaranteed.

- Revolutionary tissue sample fixation container:
- Guarantees zero exposure to formalin.
- No spills! No harmful exposure to gases!
- Streamlines the vial preparation process.
- Already a standard solution in many countries in Europe and overseas.



## THE SECRET IS THE LID



Twist off the lid



Place biopsy in the vial



Twist on & tighten the lid



Press of the thumb to realease formalin



The formalin flows to the biopsy



Biopsy is ready for SAFE transportation

ART.CODE	DESCRIPTION	UNITS / BOX
1748	Safe Biopsy container with 4% buffered formaldehyde 20ml	Box of 24pcs.
1749	Safe Biopsy container with 4% buffered formaldehyde 60ml	Box of 18pcs.